

All singing, all dancing
multi-positional level 4 care chair

For those with complex needs, suitable for home,
healthcare, nursing and residential environments.

Cabaret




BRITISH DESIGN - BRITISH MANUFACTURE

Cabaret

Multi-positional level 4 care chair Standard Features


Back rest recline angle,
with interchangeable seat
and back cushions.



Tilt in space action

Removable seat cushion
modules for maximum
pressure relief.
Takes any cushion

Memory foam surfaces
for pressure relief on all
contact surfaces



- a. Adjustable arm height 6" to 12"
- b. Adjustable seat depth 18" to 22"
- c. Adjustable seat width 16" to 21"

Angle adjustable foot plate
that is height adjustable and neatly
folds away when not required.

With removable cover.

Extending individually
controlled leg rest

Height adjustable support tubes for thoracic supports and wings

Electric operation of leg rest and tilt in space by handset

Push handle including lever for backrest angle recline

Handle for depth adjustment

Battery back up as standard (30 cycles)

Locking knob for seat depth adjustment

Braked rear castors

Swivel castors at front

Articulating back
Seven specific profiling positions to accommodate kyphosis.

Other cushioning options can be used for complex positioning requirements



Cabaret optional extras



Leg rest channel cushion



Soft head cushion



Profile head cushion



Head support cushion



Foam head rest



Adjustable wings



Lateral waterfall cushion



Profile waterfall cushion



Pelvic belts/harnesses



PrimaGel cushion



Cool Gel cushion



Air cushion



80 cycle battery back up

Other information

Size

One size chair:-

Seat width:- 16 - 21"

Overall chair width:- 28"

Fabric options

Can be made in any Crib 5 stretch fabric to enhance pressure relieving properties on all contact surfaces.

Outer vinyl/fabric to Crib 5 anti MRSA waterproof.

Brochure shows Cabaret chair in Ultima Peacock vinyl outer with Plum 4 way stretch VP cushions.

Warranty

2 year warranty

A151287 - 10/15

