

JuniorBeds

The accessory kit that makes all Accora beds suitable for children

This accessory kit is a simple, safe and cost-effective solution that is compliant with the BS EN 50637:2017 children's bed standard. It comprises of head and footboard bumpers, which reduce impact, and a safety sleeve, reducing the risk of entrapment.



Please note: side rails are an optional extra

	Kit includes	Part number
FloorBed 1 JUNKIT-0-FL1-000	Bumpers	BHMNE-0-FL4-000
	Safety sleeve	SLE-0-FL4-100
FloorBed 2 JUNKIT-0-FL2-000	Bumpers	BHMNE-0-FL4-500
	Safety sleeve	SLE-0-FL4-100
	Protection plates	PROPL-0-FL4-500
CommunityBed JUNKIT-0-CBI-000	Bumpers	BHMNE-0-FL4-000
	Safety sleeve	SLE-0-FL4-100