

DASH™

Dash 200 HD Wheelchair

The Dash 200 heavy duty chair with safety features in mind.

NOW AVAILABLE IN 26" SIZE.



CRASH TESTED

Crash tested for safely transporting wheelchair users without removing them from the chair.



STANDARD HUB BRAKES

Easy-to-use and quick-to-engage brakes for safety, particularly when travelling downhill.



QUICK RELEASE WHEELS

The wheels can be removed quickly and easily to aid in transportation and storage.

Dimensions

Total width	29"/31"/33"/35"
Total length	45.3"
Total height	39"
Seat width	20"/22"/24"/26"
Seat depth	18"
Backrest height	19.5"
Leg rest extension range	16"-20"
Front wheels	8.2"
Rear wheels	24"

Other technical data

Frame	Powder coated steel
Front wheel	PU tires
Rear wheel	Spoke wheels with PU tires, aluminium handrim and quick release
Handgrip	PVC
Seat and Backrest	Vinyl
Armpad	PVC
Footplate	Aluminum



Weight info

Chair lifting weight 25.5kg

Max user weight 200kg (20"/22"/24" models)

Max user weight 225kg (26" model)

Key features

- ✓ Crash tested
- ✓ Quick release wheels
- ✓ Optional elevating leg rest
- ✓ Handle brake
- ✓ Puncture proof tires

For enquiries:

T: 0845 146 0600

E: sales@rhealthcare.co.uk

F: 01246 453881

www.rhealthcare.co.uk

DASH

RHealthcare

It's ability that matters