

PRODUCT INSTRUCTION SHEET

REFERENCE Barkby Bariatric High Back Chair with Wings (708)
Without Wings (709)

DESCRIPTION

Brown vinyl upholstery, robust all welded metal frame, with or without winged head support and integral lumbar support. Completely adjustable (overall height, seat depth and base height). Fitted with polyurethane arm pads.

TECHNICAL SPECIFICATIONS

Seat Height 42cm (16.5") - 55cm (21")
Seat Width 61cm 24"
Seat Depth 45cm 18"
Seat to top of arm 20cm
Seat to top of back 68cm
Overall width 70cm (28")
Overall depth 65cm (26")

Weight Limit: 254kg (40st)



USAGE REQUIREMENTS

1. When adjusting the legs, ensure that there is no gap between horseshoe clip and frame. The clip spindle must be fully located right through the outside frame.
2. Check that the two wing nuts are fully tightened on the seat base and back support.
3. Do not exceed user weight of 254kg (40 stone).
4. Clean the frame with cream cleaner and upholstery with proprietary polish.