

# Ascenda<sup>®</sup>

Transfer from lying to standing. Effortlessly

Ascenda is a **'riser-bed'** designed to facilitate effortless transfers from lying to standing, enabling individuals to maintain dignity and mobility at home while reducing care package requirements.



**Accora**

# Key features



## Tilt-in-space mechanism

Enhances the individual's sense of safety when transitioning to a sleeping position.



## Rise function

Rise function to assist with sit-to-stand transfers, with integrated handles for stability.



## Sleeping position

Ascenda gently lowers to a fully horizontal position for sleeping.



## Depth-adjustable

Adjust the seat depth to 400 mm, 450 mm or 500 mm.



## Pillow clips

Clips on the bed keep pillows and cushions securely positioned during movement.



## Handset

Compact handset with intuitive interface.



Enable transfers



Enable mobility



Enable independence



Reduce care packages



Ascenda has made a big difference to our lives. I can get in and out of bed independently, improving sleep for both my wife and me. – Mr L



## Technical specification

Safe working load/Max patient weight	35 stone
Seat depths	400 mm , 450 mm or 500 mm
Seat height	400 mm - 600 mm
Ascenda mattress 25 stone capacity part number	VISMAT25-CHI-FM5-000
Ascenda mattress 35 stone capacity part number	VISMAT35-CHI-FM5-000
Part number	CHB-0-CHI-000

Contact us for more information or to arrange a demonstration.

# Accora

Charter House, Barrington Road, Orwell, Cambridge SG8 5QP T: +44 (0)1223 206 100 E: info@accora.care www.accora.care