



Buggies & Wheelbases

DYNO™



Key Features

- Folds Compactly
- 3 Wheel Configuration
- Swivelling Front Wheel
- 30° Hip Angle Adjustment
- ACTIVATOR™ Dynamic Seating
- Weighs only 13kg
- Rear Suspension
- Crash Tested

Active lifestyles need active strollers.

The Dyno provides positioning for postural support while allowing controlled dynamic movement of the child. The DYNO 3-wheel frame is built for all types of terrain. With our built-in suspension, jogging and traveling are made easy. The Dyno folds in one easy step, making life easier for families on the go.

13 - 21" Model

Weight Limit	36kg
Seat width	14" / 35 cm
Seat depth	6-14" / 15-35 cm
Back height	13-21" / 33-53 cm
Lower leg length	8-14" / 20-35 cm
Chest width length	7-14" / 18-35 cm
Folded with wheels	42x 17 x 26.5" / 107 x 43 x 67 cm
Seat to floor	85°, 90°, 95°, 100°, 105°
Hip angle adjustment	19" / 48 cm
Adjustable seat tilt	85°, 90°, 105°
Knee active range (optional)	10°
Seat to back active range (optional)	30°
Push handle weight	39-46" / 99-117 cm
Front wheel diameter	10" / 25 cm
Folded without wheels	35 x 17 x 25" / 89 x 43 x 64 cm
Rear wheel diameter	16" / 41 cm

Accessories



ACTIVITY TRAY



i2i HEADREST



COMBO SERIES HEAD SUPPORT



RESPOND™ SEATING



SUNSHADE



COMFORT PLUS HEAD SUPPORT